

AMERICA'S #1  
LIFT CHAIR BRAND



# LIFT CHAIR POSITION GUIDE

	 Lifted	 Seated	 Reading	 Napping	 Relaxing	 Sleeping	 Trendelenburg	 Zero-Gravity
<b>2-Position</b>	✓	✓	✓	✗	✗	✗	✗	✗
<b>3-Position</b>	✓	✓	✓	✓	✗	✗	✗	✗
<b>Infinite-Position</b>	✓	✓	✓	✓	✓	✓	✗	✗
<b>Zero-Gravity</b>	✓	✓	✓	✓	✓	✓	✓	✓



## 2-POSITION TV RECLINE

The 2-position TV recline is achievable on all Pride® Lift Chairs. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.



## 3-POSITION FULL RECLINE

The 3-position full recline is available on most Pride® Lift Chair models and puts you in a position ideal for napping or relaxing.



## INFINITE-POSITION SLEEP RECLINE

The infinite-position sleep recline is perfect for taking a nap – whether quick or extended – or even for a full night's sleep. Infinite-position Pride® Lift Chairs enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.



## TRENDELENBURG POSITION

The Trendelenburg position, available on the LC-525, SR-525 and LC-576, elevates one's feet higher than his or her head for increased comfort and relaxation.



## ZERO-GRAVITY POSITION

The zero-gravity position, which can be reached on infinite-position Pride® Lift Chair models, positions one's thighs and lower legs higher than his or her torso.